

MAKE	MODEL	ENT SHEATH	ENT CHANNEL	ENT SENSORY	Primary Use	Working Length	Insertion Tube OD	Channels	WC Size	Field of View	Depth of Field	Up	Down
GYRUS	800-600	22-5401	--	4-45401									
GYRUS	Duraview OL-1	22-5301	33-5201	44-5301	General	300	3,6	0	--	70°	3-50	135	135
GYRUS ACFI	N/A	22-5401	--	44-5401	General	300							
KARL STORZ	11101RP	22-6301	33-6301	44-6301	General	250	3,5	0	--	70°	N/A	180	100
KARL STORZ	11101RP-1	22-6401	33-6401	--	General	300	3,5	0	--	70°	N/A	180	100
KARL STORZ	11101RP2	22-6401	33-6401	--	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
KARL STORZ	11101SK1	22-6101	--	--	Pediatric	250	2,5	0	--	70°	N/A	180	100
KARL STORZ	11101TP	--	--	--	Therapeutic	350	3,5	0	--	70°	N/A	N/A	N/A
MACHIDA	ENT-3L	22-5301	33-5301	44-5301	General	300	3,2	0	--	80°	3-50	130	100
MACHIDA	ENT-4L	22-2401	33-2401	44-2401	General	300	4,2	0	--	80°	3-50	130	100
MACHIDA	ENT-2L	22-2101	--	--	Pediatric	N/A	2,2	0	--	N/A	N/A	N/A	N/A
MACHIDA	Videoscope	Not Tested	Not Tested	Not Tested	General	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
MACHIDA	ENT-30PIII	22-5301	33-5301	44-5301	General	300	3,2	0	--	80°			
MEDTRONIC	SharpSite	22-5401	33-5401	44-5401	General	300	3,4	0	00	80°	N/A	130	130
OLYMPUS	ENF-P3	22-4301	33-4301	44-4301	General	250	4.2*	0	--	80°	3-50	130	130
OLYMPUS	ENF-P2	22-4301	33-4301	44-4301	General	250	4.2*	0	--	80°	3-50	130	130
OLYMPUS	ENF-T3	--	--	--	Therapeutic	300	5	1	2,2	85°	2-50	120	120
OLYMPUS	ENF-P4	22-4401	33-4401	44-4401	General	300	4.2*	0	--	80°	3-50	130	130
OLYMPUS	ENF-GP	22-4401	33-4401	44-4401	Portable	300	4.2*	0	--	80°	3-50	130	130
OLYMPUS	ENF-L3	--	--	--	Therapeutic	265	4,2	0	--	85°	5-50	130	130
OLYMPUS	ENF-XP	22-4101	--	--	Pediatric	N/A	2,6	0	--	75°	N/A	130	130
OLYMPUS	ENF-V	22-4501	--	--	Video General	365	3,9	0	--	90°	N/A	130	130
OLYMPUS	ENF-VQ	22-5501	--	--									
OLYMPUS	ENF-V2	22-6401	--	--	Video General	300	3,2	0	--	90°	N/A	130	130
OLYMPUS	ENF-VT	--	--	--	Video Tx	365	4,8	1	2	90°	N/A	130	130
OPTIM	ENT 3.6	22-4401	--	--									
OPTIM	ENTity	22-4401	--	44-4401									
OPTIM	FreedomView	Not Tested	Not Tested	Not Tested	General	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
PENTAX	FNL-15RP3	--	--	--	Therapeutic	300	4,8	1	2,2	75°	N/A	130	130
PENTAX	1070STK	22-6401	33-6401	none									
PENTAX	FNL-10P2	22-5401	33-5401	44-5401	General	300	3.7*	0	--	75°	3-50	130	130
PENTAX	FNL-10S	22-5401	33-5401	44-5401	General	300	3.7*	0	--		N/A	130	130
PENTAX	FNL-13S	22-2401	33-2401	44-2401	General	300	4,1	0	--	120°		180	130
PENTAX	VNL-1330	22-2401	33-2401	44-2401	Video General	300	4,2	0	--	85°	N/A	130	130
PENTAX	VNL-1530T	--	--	--	Video Tx	300	5,3	1	2	85°	N/A	130	130
PENTAX	FNL-7RP3	22-4101	--	--	Pediatric	300	2,4	0	--	75°	N/A	130	130
PENTAX	FNL-10RP3	22-5401	33-5401	44-5401	General	300	3.7*	0	--	75°	3-50	130	130
PENTAX	VNL-1130	22-5501	--	--	Video General	300	3,8	0	--	85°	N/A	130	130
PENTAX	VNL-1170K	22-5501	--	--	Diagnostic	300	4,1	0	--	85°	N/A	130	130
PENTAX	FNL-10RBS	22-5401	33-5401	44-5401	Portable	300	3,4	0	--	75°	--	130	130
PENTAX	FNL-10RAP	--	--	--	FEESST	350	3,4	1	1,2	95°	N/A	130	130
PENTAX	FNL-13RAP	--	--	--	FEESST	350	4,1	1	1,2	100°	N/A	130	130
SCHOLLY	3.2 Sharpsite	22-5401	33-5401	44-5401									
VSI	ENT-3000	22-5401	33-5401	44-5401	General	300	3,6	0	--	80°	3-50	135	135
VSI	ENT-2000	22-5401	33-5401	44-5401	General	300	3,6	0	--	80°	3-50	135	135
VSI	ENT-1000	22-5101	--	44-5101	Pediatric	300	2,4	0	--	80°	3-50	135	135
VSI	E-F100	22-5301	33-5201	44-5301	General	300	3,6	0	--	70°	3-50	135	135
VSI	ENT-5000	22-5601	33-5601	44-5601									
VSI	ENT-4500	22-5101	--	--									
VSI	TNE-5000 2.1mm	03-5101	--	--									
VSI	TNE-5000 1.5mm	03-5102	--	--									
WELCH ALLYN	RL-150	22-5401	33-5401	44-5401	General	300	3.7*	0	--	75°	N/A	130	130
XION	EF-NS	22-6401	33-6401	--	General	300	3,6	0	--	N/A	N/A	130	100
XION	EV-N	22-5501	--	--									
XION	EV-NE	22-2401	33-2401	44-2401									

* = True Dimension, Not Marketing Number

Scope WITH channels will not have an EndoSheath